


Quick Tips:

Get the best results from your hearing aid batteries:



Due to the Covid-19 pandemic, there may be delays in getting your hearing aid batteries. Here are some tips to help extend your battery life.

1. Store batteries in a cool, dry, place.
2. Keep batteries away from coins, keys and other metals to avoid discharging the batteries.
3. Check the expiry date of batteries before use and do not use if expired.
4. Use the oldest batteries you have on hand first (those with an expiry date closest to the current date).
5. Wash your hands before handling the battery.
6. Do not remove the tab from the battery until you are ready to use the battery.
7. Wait at least one minute after removing the tab from the battery to allow the battery to activate fully before use.
8. Open the hearing aid battery door at night when not in use.
9. If you are not using your hearing aid(s) for an extended period of time remove the battery(s) entirely.
10. Remove dead batteries immediately.

Note: Hearing aid battery life is affected by a number of factors including the number of features activated, feedback, hearing aid usage, streaming and how batteries and hearing aids are cared for.